



**New!**  
**Portable Water resistant pH Meter for the Beer Industry, HI-99151.**

- **Water resistant**
- **Multi-level LCD display** (The meters display the current measurement simultaneously with the current temperature. A quickguide is also displayed to aid users in instrument operation)
- **HOLD** (The Hold button freezes readings on the display so you can easily record readings manually)
- **Automatic Temperature Compensation**
- **Automatic 1 or 2 point calibration**
- **BEPS** (BEPS (Battery Error Protection System) alerts the user in the event that low battery power could adversely affect readings)
- **Battery % displayed on startup.**
- **Heavy-duty** (The enclosure of these waterproof meters are made of tough ABS plastic and the electrode features a titanium body)

**HI 99151 is supplied with:**

- FC 214 pH electrode with internal temperature sensor
- pH 4 & pH 7 buffer sachets
- HI 700661 electrode cleaning solution (2 x 20 mL)
- batteries
- rugged carrying case
- instructions



Manufacturer's recommendations:

The Hanna HI99151 beer pH meter should be used with grab samples. Please dispose of samples after use. It is not intended to be used in tanks.

Sample temperature should not exceed 40°C / 104°

**HannaNorden AB**  
 Energigatan 15 B  
 S-434 37 KUNGSBACKA  
 SWEDEN

Phone +46 300 404018  
 Fax +46 300 14122  
 info@hannanorden.com  
 www.HannaNorden.com

